



WOLTURNUS



Basket

Rigid wheelchair
for basketball

Sebastian Magenheim with his Wolturnus Basket. Copyright: Jürgen Berg

Strong and light Basketball wheelchair

The Wolturnus Basket wheelchair is a first-class sports wheelchair. It's strong, robust and light at the same time and designed in accordance with international regulations.

The Wolturnus basketball wheelchair has been optimised in partnership with international players from e.g. German and Dutch clubs. It is designed in accordance with international regulations.

The frame is made of tempered 7020 aluminium, the strongest aluminium alloy that can be welded. The result is a first-class sports wheelchair with superior stability and low weight. It ensures that you can accelerate and switch directions faster than opponents. The frame's geometry is adapted to suit your size, strength and position on court.

Because the frame is anodised, the surface is hard, robust and maintenance-free.

The Wolturnus Basket is available as a fully welded customised version as well as a model with adjustable features. Basket A is therefore suitable for beginners as well as a team wheelchair.


Basket key features

- › Standard and customised models
- › Fully welded or with adjustability features
- › Design adapted to suit your position on court
- › Extremely easy to manoeuvre and to achieve speed
- › High stability and low weight
- › Silver anodised 7020 aluminium frame; corrosion resistant, impact resistant, scratch resistant
- › Designed in accordance with international regulations
- › Danish design and manufacture



› BASKET

Customised basketball wheelchair for experienced users. Fixed backrest and rear axle.



› REAR AXLE

You can choose between a fixed and an adjustable rear axle.



› ANTI-TIP

Anti-tipping device in both sides for maximum stability on court