



WOLTURNUS

Rex 350

Power Wheelchair



Reliable and durable power wheelchair

The Rex 350 is a very stable and slim power wheelchair. It has a top speed of max. 14 km/h and a range of 45 km.

The Rex 350 is a compact power wheelchair for children and adults. This wheelchair model has a small turning radius, it is easy to manoeuvre and very quiet. It is suitable for indoor and outdoor use. The Rex 350 has a very low seat height, which makes it easy wheel in under tables.

The Rex 350 offers a good variety of seating positions. There are up to five different electrically-operated seating positions like e.g. elevated leg supports, an adjustable back angle, a seat lift as well as a seat tilt.

As an optional extra, there is a third-party steering available for situations, where the user is not able to operate the wheelchair independently.

The seat modul of the Rex 350 can be customised based on each individual user's requirements and needs. The Rex 350 is crash-tested and can be used as an in-car seat.

It is possible to equip the Rex 350 with the Wolturnus Wing Back ILSA.

Key features

- › Compact, reliable and durable power wheelchair
- › Easy to manoeuvre and very quiet
- › Impact-resistant anodised topframe
- › Total width only 58 cm
- › Top speed of max. 14 km/h
- › Individually fitted seat unit
- › Possibility for adjustable strap seat & strap back
- › Up to five electrically-operated seating-positions
- › Crash-tested
- › Can be used as in-car seat
- › Weight from 84 kg
- › User weight max. 140 kg



› LIFT & TILT

The chair is available with electrically-operated seat lift and seat tilt



› PUNCTURE-PROOF WHEELS

As a standard, Rex 350 comes with 12,5" puncture-proof rear wheels and 7" puncture-proof castors.



› ACCESSORIES

There are a lot of accessories available to choose from. E.g. headrest, armrests, Wing Back ILSA etc.