



SYST'AM® P913L / OCCIPITAL SYSTEM

PREVENTION AND POSITIONING OF THE OCCIPITAL AREA

MATERIALS

80
kg/m³



Visco foam



Removable
POLYMAILLE®
integral cover

NF EN ISO 597 - 1 & 2

Foam maintenance:



Cleaning POLYMAILLE® cover:



INDICATIONS

- Prevention and positioning of the occipital area.



VISCOELASTIC FOAM WITH SHAPE MEMORY

- Defined body moulding and increased body surface area in contact with the foam.
- More effective distribution of pressure: reduction of transcutaneous pressure on highrisk zones, facilitated blood flow.
- Increased comfort and more effective maintenance of patient position.

AVAILABLE VERSIONS



SYST'AM® P913L /
OCCIPITAL SYSTEM

FEATURES OF THE COVERS

FABRIC COATED WITH BI-STRETCH POLYURETHANE

- Reduces friction and shear effects.
- Supple and soft to the touch (comfortable).
- Favours the exchange of gases (steam, sweat):
 - fights against maceration.
- Impermeable material:
 - better hygiene,
 - longer support system lifespan.
- Washable at 90°C, can be decontaminated using cold sprays.
- Treated to resist fire.
- In multi-patient use, it is preferable to buy one new cover per patient.



DESIGNATION	ITEM CODE	SIZE L x W x H (cm / inches)
Occipital system	P913L1HW	50 x 43 x 9,5 cm / 19,7 x 16,9 x 3,7"

