



WOLTURNUS

# Rugby

Specialized  
sports wheelchair



# Rugby Attack & Defence

## For high speed and hard tackles

The Rugby Attack and the Rugby Defence are specialized rugby wheelchairs. They have been optimized in accordance with the most up-to-date input from Danish, Dutch and German top players.

The raw frame of heat-treated 7020 aluminium has the best strength-to-weight ratio on the market. It transfers 100% of the player's strength and capacity to the playing field.

The Wolturnus rugby wheelchairs are fast and stable. They make it easy to accelerate and switch directions. The fully welded frames are extremely strong and able to withstand repeated and hard tackling.

Both models are equipped with strong anti-tippers to secure safe manoeuvring in the heat of the moment.

The rugby wheelchairs are designed and manufactured in Denmark. They are customised to suit the individual player and the geometry of the frame is adapted to suit your strength, size and playing position.


### Key features

- › Customized to suit your requirements and wishes
- › Two basic models that are adapted to suit your handicap and playing position
- › Extremely easy to manoeuvre and to achieve high speed
- › High stability
- › Fully-welded frame in raw, heat-treated, lightweight aluminium 7020 T6
- › Designed in accordance with international regulations



#### › RUGBY ATTACK

Reinforced side bumpers to withstand hard tackles.



#### › RUGBY DÉFENSE

The extended front bumper is designed to stand your ground and block your opponent's wheelchair.



#### › ALUMINIUM 7020 T6

The raw frame of heat-treated 7020 aluminium has the best strength-to-weight ratio on the market.